



Cleaning and Sanitising Your Brewing Equipment

Cleansers and sanitizers are not the same: **cleansers** leave your equipment bright and clean; **sanitizers** keep bacteria and other spoilage organisms under control. Before you sanitize, everything must be clean. Rinse bottles, primary fermenters and carboys as soon as they are empty to avoid scraping dried, caked-on residue off them later.

Ale Supply LLC

Product	Description	Cleaning	Sanitizing	Cautions
TRISODIUM PHOSPHATE (Saniton)	Trisodium Phosphate (TSP) is an unscented detergent. It helps take a lot of the work out of removing stains and stubborn grime from equipment. TSP is especially useful for stainless steel, which a chlorinated cleaner such as Diversol (discussed below) could corrode.	<ol style="list-style-type: none"> 1. Dissolve 2 g per litre of warm water (2 tsp. per gal). 2. Soak equipment for 20 minutes; scrub any stubborn stains. 3. Rinse well with hot water. 	TSP does not sanitize. You must use it in conjunction with a proper sanitizer.	Avoid contact with strong acids and pro- longed contact with aluminum, tin, lead, and zinc since this will produce hydrogen gas.
DIVERSOL (Sani-Brew)	Diversol is extremely effective both as a cleaner and as a sanitizer. Available under a variety of trade names, you'll recognize Diversol because it is a chlorinated pink powder.	<ol style="list-style-type: none"> 1. Dissolve 3.5 g per litre of cold water. (5 tsp. per gal). 2. Soak stained equipment for up to 48 hours. 3. Rinse thoroughly with hot water. 	<i>Sanitizing most equipment:</i> <ol style="list-style-type: none"> 1. Fill primary fermenter with Diversol solution (3.5 g per litre of cold water [5 tsp. per gal.]). Fully immerse everything you'll need: hoses, sp oons, bungs, hydrometers, thermometers, airlocks, and primary lid. 2. Soak equipment for at least 20 	Corrodes stainless steel. Can bleach clothing. Do not mix with acids, amines, or ammonia. Such a mixture produces dangerous gases.

TI

Quality Wine and Ale Supply ~ www.HomeBrewIt.com

“Elkhart County’s First and Finest Winemaking and Homebrew Specialty Store”

Store: 108 S. Elkhart Ave., Elkhart, IN 46516 ~ Mail: 530 E. Lexington Ave., Suite 115, Elkhart, IN 46516 USA

Website: www.HomeBrewIt.com ~ 24/7 Email: info@HomeBrewIt.com ~ 24/7 Phone/Fax/Voice: +1 (574) 295-9975

Prices and policies are subject to change without notice. All rights reserved © Quality Wine and Ale Supply LLC

The contents of this document is the property of the copyrighted owners and is provided ‘as is’ without guarantee of accuracy.



			<p>minutes. 3. Remove lid, rinse carefully (but thoroughly) with hot water, and turn it upside-down on the counter. Rinse equipment, placing items inside sanitized lid. (Fig. 1)</p> <p><i>Sanitizing a carboy:</i> 1. With bung in carboy, slosh 9 litres (2 gal.) Diversol solution all around inside. (Fig. 2) 2. Repeat twice, with 5 minutes between repetitions. Rinse after 20 minutes.</p>	
<p>BLEACH (Sodium Hypochlorite)</p>	<p>Unscented household bleach can be used as a sanitizer and cleaner. For sanitizing, it is just as effective as Diversol. For cleaning, however, either Diversol or TSP would be preferable.</p>	<p>Mix 2.5 ml per litre of cold water (1 tbsp. per gal.) and follow the cleaning instructions for Diversol (above). When rinsing, remove all traces of chlorine smell.</p>	<p>Using a solution of 5 ml per litre of cold water (1 tbsp. per gal.) follow the cleaning instructions for Diversol (above.)</p>	<p>See above. Don't use scented bleach. Those "spring fresh" and "lemony" perfumes will sink into equipment permanently. All subsequent batches will smell and taste like perfume.</p>
<p>IODOPHOR</p>	<p>If you keep your equipment clean between brews, Iodophor is a great sanitizer. (Your equipment will not be clean if you simply rinse after use. You</p>	<p>Iodophor does not clean. Use a cleanser such as TSP before sanitizing with Iodophor.</p>	<p>1. Mix 1 ml per litre of cold water. You must use a syringe for accurate measurement. 2. Use a spray bottle to coat equipment with solution.</p>	<p>When in solution, Iodophor rapidly breaks down. Every time you use it, you must make a fresh solution. Not effective</p>

TT

Quality Wine and Ale Supply ~ www.HomeBrewIt.com

“Elkhart County’s First and Finest Winemaking and Homebrew Specialty Store”

Store: 108 S. Elkhart Ave., Elkhart, IN 46516 ~ Mail: 530 E. Lexington Ave., Suite 115, Elkhart, IN 46516 USA

Website: www.HomeBrewIt.com ~ 24/7 Email: info@HomeBrewIt.com ~ 24/7 Phone/Fax/Voice: +1 (574) 295-9975

Prices and policies are subject to change without notice. All rights reserved © Quality Wine and Ale Supply LLC

The contents of this document is the property of the copyrighted owners and is provided ‘as is’ without guarantee of accuracy.



must use a cleaner, such as TSP.) Iodophor is excellent for stainless steel- which Diversol will corrode- and all other beermaking equipment. Mixed properly and used in a spray bottle, Iodophor is almost magical in its convenience: there's no waiting and no rinsing.

3. Shake off excess or allow equipment to drain. Do not rinse.
 4. Use the equipment. **It's that easy!**

unless mixed accurately. Use a syringe. Will not clean dirty or stained equipment. Hazardous for those with iodine allergies. Concentrated solution stains skin and fabric.

FINAL NOTES

If you keep your equipment looking, smelling and feeling clean, even when it's not in use, you'll have fewer problems keeping everything sanitized.

Store your clean primary fermenters and carboys closed with a little chlorine solution inside (one capful unscented household bleach per litre cold water). Nothing will grow in them, so sanitizing will be a snap.

If you get tired of scrubbing and sanitizing, remember: lapses in sanitation are responsible for 90% of all homebrewing failures. You can avoid them.

This document is provided by Quality Wine

Quality Wine and Ale Supply ~ www.HomeBrewIt.com

“Elkhart County’s First and Finest Winemaking and Homebrew Specialty Store”

Store: 108 S. Elkhart Ave., Elkhart, IN 46516 ~ Mail: 530 E. Lexington Ave., Suite 115, Elkhart, IN 46516 USA

Website: www.HomeBrewIt.com ~ 24/7 Email: info@HomeBrewIt.com ~ 24/7 Phone/Fax/Voice: +1 (574) 295-9975

Prices and policies are subject to change without notice. All rights reserved © Quality Wine and Ale Supply LLC

The contents of this document is the property of the copyrighted owners and is provided ‘as is’ without guarantee of accuracy.